

Porous Quartz Filter Discs



The key to success when working with quartz glass filters is selection of the correct porosity. Glass filters vary in porosity according to grades that are labeled 0 to 5. The table below shows the porosity ranges and the main fields of application for each grade. The pore size indicated always refers to the largest pore in the disc. This also indicates the diameter of the smallest particle retained during filtration. Porosities are determined using the Bechhold bubble-pressure method which has often been described in literature.



Porosity	Nominal max. pore size (μm)	Fields of applications, examples
00	250 - 550	Liquid and gas distribution
0	160 - 250	Gas distribution Gas distribution in liquids at low pressure Filtration of very coarse precipitates
1	100 - 160	Coarse filtration, Filtration of coarse precipitates, gas distribution in liquids Liquid distribution, coarse gas filtration Extraction apparatus for coarse grain materials Loose filter layer substrates for gelatinous precipitates
2	40 - 100	Preparatory fine filtration Preparatory work with crystalline precipitates. Mercury filtration.
3	16 - 40	Analytical filtration. Analytical work with medium-fine precipitates. Preparatory work with fine precipitates. Filtration in cellulose chemistry, fine gas filtration Extraction apparatus for fine-grained materials
4	10 - 16	Analytical fine filtration Analytical work with very fine precipitates (e. g. BaSO_4 , Cu_2O) Preparative work with precipitates of appropriate fineness. Non-return and stop valves for mercury

Flow rates

Water and air flow rate through filter discs of various porosities as a function of pressure differential. For disc diameter 30 mm.



Ultrafine Filtration

Ultrafine filtration is one of the most important methods for treatment of biological solutions without using high temperatures which often lead to changes in or decomposition of the active ingredients in the solution. For liquid filtration, sintered glass filter funnels of standard design are used. Dilute suspension of bacteria (15 000 to 90 000 per millilitre) can be sterilely filtered using porosity 3 filters. The pores are already so narrow that all bacteria in dilute suspensions adhere to the pore walls. However, a bacterium-free filtrate could not be obtained when filtering dense suspensions through them. In the case of dense suspensions, bacteria can still pass through once the pore walls become saturated. Real straining is only achieved with a maximum pore size of 2 μm and below; i.e. it is only here that the pores are smaller than the bacteria to be separated.

For **bacterium-free filtration of gases**, e.g. in ventilation of fungal and bacterial cultures, pipeline filters are used. In this case, porosity 3 is adequate provided that the space in front of the dry filter disc on the air inlet side is stuffed evenly and loosely with cotton wool.

Standard filter discs

Other dimensions, porosities and tighter tolerances upon request

Round

OD	OD Tol	Porosity	Thickness	Thickness Tol
5	±0,5	1-4	2-4	±0,5
10	±0,5	00/0-4	2-4	±0,5
15	±0,5	00/0-4	4-8	±0,5
20	±0,5	00/0-4	4-8	±1,0
25	±0,5	00/0-4	4-8	±1,0
30	±0,5	00/0-4	4-8	±1,0
35	±0,5	00/0-4	4-8	±1,0
40	±0,5	00/0-4	4-8	±1,0
45	±0,5	00/0-4	4-8	±1,0
50	±0,5	00/0-4	5-8	±1,0
55	±0,5	00/0-4	5-8	±1,0
60	±0,5	00/0-4	5-8	±1,0
70	±0,5	00/0-4	5-8	±1,0
80	±0,5	00/0-4	5-8	±1,0
90	±0,5	00/0-4	5-8	±1,0
100	±1,0	00/0-4	8-30	±2,0
110	±1,0	00/0-4	8-30	±2,0
120	±1,0	00/0-4	8-30	±2,0
130	±1,0	00/0-4	8-30	±2,0
140	±1,0	00/0-4	8-30	±2,0
150	±1,0	00/0-3	8-30	±2,0
200	±1,0	00/0-3	8-30	±2,0
250	±1,0	00/0-2	8-30	±2,0
300	±1,0	00/0-2	8-30	±2,0

Rectangular

up to length	up to width	L, W Tol	Thickness	Thickness Tol	Porosity
400	355	±1,0	20	±2,0	00/0-2

Germany

Heraeus Quarzglas GmbH & Co. KG

Base Materials

Kleinostheim, Germany

Phone +49 (6181) 35-7444

sales.hqs.basematerials.de

@heraeus.com

USA

Heraeus Quartz North America LLC

Base Materials Division

Chandler, AZ 85226, USA

Phone (+1) 512 989 05 03

sales.hqs.basematerials.us

@heraeus.com

UK

Heraeus Quartz UK Ltd.

Neptune Road

Wallsend, Tyne & Wear

NE28 6DD, United Kingdom

Phone +44 (0)191.2598411

sales.hqs.tsc.uk@heraeus.com

www.heraeus-quarzglas.de