

ADDITIONAL INFORMATION

Caution:

Federal law restricts this device to sale, distribution, and use by or on the order of a physician.

Warnings:

If the **PALAMIX**[®] cement gun has been sterilised with heat or steam, it must cool down before reuse.

The **PALAMIX**[®] cement gun should only be used by trained personnel. The risk of hand injury or damage to protective gloves is possible due to improper handling. There is also a risk of damage to respective mixing cartridge if personnel are not trained correctly.

Precautions:

It is recommended that the gun grip bearings are oiled after each cleaning cycle.

Check the functionality of the **PALAMIX**[®] cement gun prior to each operation. If damaged, contact your local Heraeus Medical representative. The user must ensure the correct use of the **PALAMIX**[®] cement gun.